


**BINGO-Evolution tilbehør**

Mål 2017


**HOGGI<sup>®</sup>**

## Mål kurv modell 3 kg




Ca. mål i (cm)		A	B	C
3231-9401E	Str 1	40	24	17
3233-9401E	Str 2	45	24	17

## Mål kurv modell 8 kg




Ca. mål i (cm)		A	B	C	D
3231-9402E	Str 1	37	32	24	15
3233-9402E	Str 2	42	32	24	15

## Abduksjonskloss




Ca. mål i (cm)		A	B	C
	Str 1	9	6	9
	Str 2	11	8	10
	Str 3	12,5	9	10

## Bryst-skuldersele




Ca. mål i (cm)		A	B	C	D	E	F
3201-5600	Str 1	18	9	8,5	18,5	27	11
3201-5700	Str 2	23,5	12	10	21,5	35	15
3203-5600	Str 3	24	12	11	28	41	17

## Vest




Ca. mål i (cm)		A	B	C	D	E	F
3201-5800	Str 0	21	25	11	8,5	9	4
3201-5400	Str 1	32	31	11	11	10	5
3201-5500	Str 2	37	39	14	13	15	5
3203-5400	Str 3	42	45	16	13	17	6

## Hodestøtte, stor



Ca. mål i (cm)		A	B	C	D	E
3201-7600 3201-7601	Str 1	17	10	5	8 - 17	12 - 21
3203-7600 3203-7601	Str 2	19	11	5	13 - 22	17 - 26

## Hodestøtte, liten



Ca. mål i (cm)		A	B	C	D
3231-7600	Str 1	16	6	14 - 16	20 - 22
3233-7600	Str 2	18	7,5	18 - 21	22 - 25




## Sidestøtter




Ca. mål i (cm)		(A)	(B)
3231-7700	Str 1	8	10
3233-7700	Str 2	10	12
3234-7700	Str 3	12	13
Avstand (C)	Sete str 1	Sete str 2. / Sete str 2XL	
		12 - 22	17 - 27

## Bord



Ca. mål i (cm)		A	B	C	D
323x-8200 323x-8202	Plexi	52	47	29	16
323x-8201 323x-8203	Grå plast	52	47	29	16

## Bord høyderegulerbar



Ca. mål i (cm)		A	B
3231-820X	Str 1	16,5	13,5
3233-820X	Str 2	18	15
3231-820X (Høyde min-max)	Str 1	max. 21 min. 17,5	max. 18 min. 14,5
3233-820X (Høyde min-max))	Str 2	max. 22 min. 18,5	max. 19 min. 15,5


# HOGGI®

## Varmepose



Ca. mål i (cm)		A	B	C	D	E
3201-8300	Str 1	38	53	60	46	19
3201-8400	Str 2	43	60	74	60	20
3203-8300	Str 3	48	67	90	76	21

## Mål basisplate-spesialunderstell for ekstra utrustning



Ca. mål i (cm)		A	B	C	D
3201-4100	Str 1	27	18	38,5	29
3203-4100	Str 2	32	18	43,5	29

max last 15 kg

**Hersteller/Manufacturer:**

**HOGGI GmbH**

Eulerstraße 27

D - 56235 Ransbach-Baumbach

Fon: (+49) 26 23 / 92 499-0 · Fax: (+49) 26 23 / 92 499-99

e-mail: [info@hoggi.de](mailto:info@hoggi.de) · Internet: <http://www.hoggi.de>