

Vedr.: Nyt hjul til Mini Crosser X	PUP nr.:	21-007
	Udarbejdet af:	LPA
	Dato:	17-03-2021
	Intern <input type="checkbox"/>	Ekstern <input checked="" type="checkbox"/>

Vi lancerer et nyt hjul med grovere dæk, til Mini Crosser X-model, så den står endnu bedre fast på underlaget.

Varenr.:

1503-1016 HJUL, MODEL-X, 2,50-3,70-8" GROV, KOMPLET MED FÆLG OG DÆK

1503-1396 HJUL, MODEL-X, 2,50-3,70-8" GROV, KOMPL. PUNKTERFRI PUR

1503-1400 HJUL, MODEL X, 2,50/3,70-8", GROV, KOMPLET PIGGET

1503-1402 HJUL, MODEL-X, 2,50-3,70-8" GROV, KOMP MED FÆLG, DÆK, PANSERRING

Reserve dele:

1503-1397 DÆK MODEL X, 2,50/3,70-8", GROV, PIGGET DÆK

1503-1398 DÆK MODEL X, 2,50/3,70-8", GROV

Slange og fælg er det samme som standard hjulene.



Grov komplet med fælg og dæk



Grov pigget dæk



Klaus Smedegaard
Administrerende direktør



Regarding: New wheel for Mini Crosser X	PUP No.:	21-007
	Prepared by:	LPA
	Date:	17-03-2021
	Internal <input type="checkbox"/>	External <input checked="" type="checkbox"/>

We are launching a new wheel with coarser tires, for the Mini Crosser X model, so that it sticks even better to the ground.

Art.No:

1503-1016 Wheel, model-X, 2,50-3,70-8" coars, complete with rim and tyre

1503-1396 Wheel, model -X, 2,50-3,70-8" coars, complete. Puncture free, PUR

1503-1400 Wheel, model X, 2,50/3,70-8", coars, complete with spikes

1503-1402 Wheel, model -X, 2,50-3,70-8" coars, complete with rim, tyre amd, armor ring

Reservele:

1503-1397 Tyre model X, 2,50/3,70-8", coars, with spikes

1503-1398 Tyre model X, 2,50/3,70-8", coars

Tube and rim are equal to the standard wheels.



Coars wheel w/rim and tyre



Coars tyre with spikes



Klaus Smedegaard
General Manager

