

Vedr.:  Information om gaspotentiometer M & X model fra 2011	PUP nr.:	25-002
	Udarbejdet af:	LPA
	Dato:	11-02-2025
	Intern <input checked="" type="checkbox"/>	Ekstern <input checked="" type="checkbox"/>

Ved udskiftning af gaspotentiometer på M & X modeller fra 2011 og frem, indeholder sættet 2 sæt skruer til montage af wig wag.

Til X1/M1 bruges størrelse M4x16 (8031-1292)

Til X2/M2 bruges størrelse M4-12 (8031-1515)

I sættet medfølger også denne info-seddel.



Lige nu er M4x12 blanke, men senere bliver de udskiftet med en sort model.

Alle reservedele der udsendes fra nu, vil indeholde de to sæt skruer.



**Klaus Smedegaard**  
Administrerende direktør



Regarding:  Information about throttle potentiometer M & X model from 2011	PUP No.:	
	Prepared by:	LPA
	Date:	11-02-2025
	Internal <input checked="" type="checkbox"/>	External <input checked="" type="checkbox"/>

When replacing the throttle potentiometer on M & X models from 2011 onwards, the kit contains 2 sets of screws for mounting the wig wag.

For X1/M1 use size M4x16 (8031-1292)

For X2/M2 use size M4-12 (8031-1515)

This information sheet is also included in the kit.



Right now the M4x12 are blank, but later they will be replaced with a black model.

All spare parts sent from now on will contain the two sets of screws.



**Klaus Smedegaard**  
General Manager

